

PECF EPFL

Taq polymerase isolation from *E. coli*

Solutions to prepare:

Bacterial growth media

LB

Turbo/2XY or Terrific broth

Ampicillin 1000x stock 100mg/ml

IPTG 1000x stock 500mM stock

Lysozyme

buffer A: 50mM Tris pH8, 50mM Glucose, 1mM EDTA

buffer B: 10mM Tris pH8, 50mM KCl, 1mM EDTA, 0.5% Tween 20, 0.5% Nonidet P-40 or IGEPAL

10x PCR buffer: 500mM Tris pH 9.2, 17.5mM MgCl<sub>2</sub>, 150mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>

- Culture an overnight stock with glycerol (10ml LB + ampicillin 100 $\mu$ g/ml)
  - Inoculate 200ml with Turbo broth (or a rich medium 2XY or Terrific broth + ampicillin 100 $\mu$ g/ml) with the 2ml o/n culture at 37°C until the OD 600= 0.6
  - Induce the Taq expression with IPTG 0.5mM and incubate 16h at 250rpm 30°C (37°C)
  - Centrifuge the bacteria for 15' in 50ml tubes
  - Resuspend in 11.2ml (4 tubes x 2.8ml/tube) buffer A
  - Transfer to 15ml Falcon tube
  - Prepare 2ml of buffer A, which contains 64mg Lysozyme
  - Add 0.4ml in each tube (Lysozyme 4mg/ml final)
  - Mix by inversion
  - Incubate 15' RT
  - Add 12.8ml (3.2ml/tube) buffer B
  - Mix by inversion
  - Incubate 1h at 75°C in a water bath and agitate time to time
  - Centrifuge at max speed for 20' 4°C
  - Keep 5 min on ice
  - Take supernatant, and add glycerol (50/50 with water)
  - Freeze at -80°C
  - Without Glycérol, stable for 1month à 4°C
- 
- From the glycerol stock, titrate against commercially available Taq (range of 0.5 ~ 4 $\mu$ l per PCR reaction)